



# D-Check

## Bio-Acoustic Simulator

The D-Check Bio-Acoustic Simulator performs dependable daily calibration checks of the SMART TONE Automatic Audiometer, and can be used for OSHA, MSHA, or DOT hearing programs.

Because the D-Check is powered by the SMART TONE audiometer, there is no need to worry about battery life. This advanced design simplifies the daily check and offers greater reliability and repeatability.

### Features:

- Quick and accurate daily calibration check.
- Uses no batteries—power supplied by audiometer.
- Color-coded ear cups reduce confusion and errors.
- Table-top or test booth wall-mount.

### Supplied Accessories:

- Patch cable
- Wall mounting bracket

### D-Check Specifications

Indicators:	Power, left response, right response
Test Frequencies:	500 through 8 kHz
Response Threshold:	65 dB HL nominal
Power:	supplied by SMART TONE audiometer
Dimensions:	7 in. x 5 in. x 4 in.
Weight:	1 Lb.



**SMART**  
DIAGNOSTIC DEVICES

20100614